

FAQs

How does NeedleBay help me remember my insulin?

Each used needle goes back into its own bay, so a quick look through the see-through lids shows which of today's doses are done.

How does NeedleBay help with handling pen needles?

Each needle is held in its bay while you fit and remove it, so you never touch it. Used needles stay locked away until they go into your sharps bin.

Which pen needles fit a NeedleBay?

NeedleBay works with the pen needles you're already prescribed, including BD Micro-Fine, NovoFine, mylife and GlucoRx.

Where can I buy NeedleBay?

Online at needlebay.com, with UK delivery for £2.95.

How long will NeedleBay last?

NeedleBay is made in the UK from tough, durable materials, designed for daily use, year after year.

How do I clean NeedleBay?

In warm soapy water or the dishwasher.

Does NeedleBay need batteries or an app?

No. NeedleBay is purely mechanical: no batteries, no charging and no app.

FEATURES



Greater control of your medication



See your doses at a glance



No need to touch a needle



Independent living



Single-use needles



Safe needle storage



Designed to last



Perfect for on the go



All ages

NeedleBay[®]
Medical Technology

NeedleBay Systems Limited
24-28 St Leonards Road
Windsor, Berkshire SL4 3BB
United Kingdom
0345 835 8130
info@needlebay.com



FEATURES AND BENEFITS

Needle Bay[®]
Medical Technology

Daily Needle Control

A simple, mechanical way to store your insulin pen needles and see at a glance which of today's doses are done.



NeedleBay has made many of my patient's lives easier and supports safe & hygienic injection practice. For one patient in particular it has meant that district nurses no longer have to visit 2-3 times daily to administer insulin.

Karen Pudge
NHS Diabetes Specialist Nurse



NeedleBay[®]

Medical Technology

NeedleBay helps people who inject insulin keep track of their day, with no apps, batteries or charging.

NeedleBay is a small, purely mechanical case with four bays. Each bay holds one pen needle in its outer cap, with the seal intact. Your pen screws onto a fresh needle in its bay, and after your injection the used needle goes back into the same bay.


A quick look through the see-through lids shows which of today's doses are done. A used needle has no seal, so you can tell by touch too.

NeedleBay is designed in close consultation with healthcare professionals and people living with diabetes, and made in the UK.

Up to 88%

In a trial at Cardiff University School of Healthcare Sciences with over 250 participants, NeedleBay showed up to an 88% improvement in self-management of diabetes medication.

HOW IT WORKS




1 Open the lids and flip down the click-lock.



2 Pop in your pen needles, each in its outer cap with the seal on, and close the click-lock to hold them in place.



3 When you're ready for your insulin, open the lid of a fresh needle and peel off the seal.




4 Screw your pen onto the needle and lift it out. Your pen is ready to use.



5 After your injection, screw the pen back into its bay and unscrew. The used needle stays behind in its bay.



6 Look through the see-through lids to see which of today's doses are done.



7 When the bays are full, open the lids over your sharps bin and release the click-lock. The used needles drop straight out without being touched.