## **FAQs**

How does NeedleBay help to remind me to take my insulin injection? With the pen needles clamped into the NeedleBay Pod you can easily see whether a pen needle has or has not been used.

How is NeedleBay safer for storing and using pen needles? Once clamped in the pod, there's no need to handle the pen needle during attachment or removal of an injection pen. Used pen needles are securely stored until disposal.

Can NeedleBay hold any make of pen needle? Yes, most of the leading brands can be used with NeedleBay.

Where can I buy NeedleBay? NeedleBay products can be purchased from pharmacies or direct from needlebay.com where you will also find a list of suppliers and stockists.

How long will NeedleBay last? NeedleBay has been manufactured from high quality materials and has been designed for continuous use over several years.

**How can I clean NeedleBay?** NeedleBay can be cleaned in warm soapy water or placed in a dishwasher.

**Is NeedleBay recyclable?** Yes, manufactured using 100% renewable energy, NeedleBay is committed to sustainability, with its pods being 100% recyclable.

## **FEATURES**



**Greater Control of Medication** 



Reminder System & Dose



**Proof Needle Handling** 



Independent Living



Single use Needles



Needle Storage



Cost Effective



Perfect for on-the-go



All Ages





needlebay.com needlebay.com needlebay.com



NeedleBay has made many of my patient's lives easier and supports safe & hygienic injection practice. For one patient in particular it has meant that district nurses no longer have to visit 2-3 times daily to administer insulin.

Karen Pudge, NHS Diabetes Specialist Nurse

NHS



Brilliant, quick and easy.
NeedleBay is the best!
Les Cobbett, 82



## Needle Bay Medical Technology

NeedleBay supports individuals to control their diabetes, empowering them to live life to the fullest.

NeedleBay offers a ground breaking solution for effortless insulin administration, making needle handling and storage safe and secure for use by those with diabetes.


NeedleBay stands alone in its innovation, transforming lives with its straightforward yet reliable system for administering insulin at home or on the go. Once loaded, NeedleBay ensures precise insulin doses with each use and safely contains the needle, eliminating direct contact for the user. Its compact design provides clear proof of insulin usage, essential for effective blood glucose management.

Manufactured using 100% renewable energy, NeedleBay is committed to sustainability, with its pods being 100% recyclable.

88%

FACT: NeedleBay users said that NeedleBay improved their medication and injection management by 88%.

## **HOW IT WORKS**



Open the Pod lids and flip down the pen needle click-lock.



Pop in your pen needles and close the click-lock to ensure your needles are held securely in the NeedleBay pod.



When you're ready to take your insulin, flip open the lid containing the new pen needle and peel of the seal.



Screw in your insulin pen and remove. Your pen is now ready to use.



Once used place the pen back into the Pod, and unscrew and remove the pen leaving your used needle safely in the pod.



The clear lids enable you to easily see the injections you've already had today.



When you are ready to dispose of the needles, open the lids and release the clicklock, the used pen needles can now drop safely into your sharps bin.